

Ramps



2 m model



Telescopic aluminum ramps with non-slip surface for maximum safety.






3 m model



Customized Vermeiren products

Vermeiren allows you to create your own customized ramp system.
Kindly contact our sales department.



Model	Length (mm)	Width (mm)	Height (mm)	Weight (kg)
691	1060 - 1900	250	50	7 x 2
691	1660 - 3130	250	50	10 x 2

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

