

Sedna PREMIUM

Always mobile



Sedna scooters are perfect to carry along as they are foldable in a few seconds. They are designed with a removable lithium battery with anti-theft security and inflatable tires. Equipped with a swivel seat and collapsible armrests



VERMEIREN
we care for you

Transport

The **Sedna Premium** scooter is the perfect transport scooter. Lightweight, compact in size and foldable under a minute, super easy to carry-on. Ideal for short distances and handy to negotiate sharp turns in small rooms, restaurants or shops.

Colour



Optional shopping bag under the seat



Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Sedna Premium** complies with the highest norms, laid down by certified test institutes in most European countries.



Technical data

Total width	530 mm
Total height	870 mm
Total length	960 mm
Seat height	540 mm
Turning radius	820 mm
Battery	24V - 10,5 Ah

Driving range	14 km*
Motors	270 W
Speed	6 km/h*
Weight	31,6 kg
Maximum load	115 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenpaleis 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com

For more information, contact your local dealer

