

175bis

TOILETCHAIR



total width	585 mm
seat depth	440 mm
height armrests	230 mm
seat height	460 mm
total height	805 mm
total length	550 mm
height back	355 mm
weight	9,15 kg
frame color	C26 - dark green
upholstery	black

Deviations in sizes and weight (1,5 mm, 1,5 kg, 1,5°) are possible. Subject to technical alterations.



N.V. Vermeiren N.V. • Vermeirenplein 1-15 • B-2920 Kalmthout • Belgium

Tel.: +32 (0)3 620 20 20 • Fax: +32 (0)3 666 48 94

info@vermeiren.be • www.vermeiren.com