

139

Commode chair



Model **139** is a high-quality wheeled commode chair combining easy-to-clean upholstery with a reliable and sturdy frame. The one-handed bucket action makes care more hygienic and efficient for the attendant. The swing-away armrests facilitate easy transfer.



VERMEIREN
we care for you



Comfort armpads



Swing away and removable footrests with retractable footplates



Wheels with brakes for safety



One-handed bucket operation

Colours

Frame

blue

Upholstery

black

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **139** comply with the highest norms, laid down by certified test institutes in most European countries.



The swing-away armrests facilitate easy transfer.

Technical data

■ Total width	540 mm
■ Total height	950 mm
■ Total length	820 mm
■ Seat height	528 mm
■ Seat width	420 mm
■ Height armrests	225 mm

■ Height backrest	390 mm
■ Seat depth	420 mm
■ Weight	13 kg
■ Maximum user weight	120 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer