

265 - 266 - 267

WALKING FRAMES

Model 265 demountable

Tubular steel frame painted in blue.
Padded armpit supports which are adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.



C27
blue




Model 266 foldable

Tubular steel frame painted in blue. Padded armpit supports which are adjustable in height. 4 wheels of 75 mm of which 2 with brake. Foldable for easy transport and storage.



Model 267 demountable

Same features as the mod.265 with the addition of the adjustment of the base height up to 150 mm.



265	1170-1470	760	520	580	min. 300	13,7	120
266	1170-1470	760	520	580	min. 275	12,5	120
267	1115-1560	760	465-610	580	min. 300	14	120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

