

# Lift chairs

Ontario 1 - Ontario 2  
Montreal

Lift chairs are recliners used to assist people in standing up. Complete with a remote control, they allow you to take a nap, watch tv, dine or use that extra boost from the lift chair to help you to get up. Available in different upholstery colors to find the perfect match with your interior.



**VERMEIREN**  
we care for you



## Ontario 1 - Ontario 2

The Ontario liftchair offers you optimal comfort in your own house. With a simple push of a button you can relax, watch tv or have a little nap in your Ontario liftchair.

The **2 motors version (Ontario 2)** allows you independent movement of the backrest and the legrest. You can adjust the legrest without the backrest and vice versa and you can flatten the liftchair completely (180°).

### Colors of the Ontario 1 / Ontario 2



*\*only for the model Ontario 1*



## Montreal

Make your choice from different colours and textures for an optimal match with your existing interior. The Montreal can be placed very close to the wall (100 mm) and still recline completely. The footrest inclines first making it easy to watch television without having to lay back completely.

### Colors of the Montreal



## Technical data Ontario 1 / Ontario 2 / Montreal

Total width	800 / 845 / 850 mm
Total height	1160 / 1145 / 1140 mm
Total length	840 / 960 / 830 mm
Seat width	465 / 480 / 485 mm
Height armrests	150 / 140 / 150 mm
Seat height	500 / 500 / 500 mm

Seat depth	450 / 530 / 510 mm
Height backrest	695 / 670 / 700 mm
Distance from the wall	450 / 450 / 100 mm
Weight	49,6 / 53,3 / 53,98 kg
Maximum user weight	150 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you