


MAT X1 X2 X3

Alternating Pressure Relief Mattresses

Ideal for:

- **Bedsore prevention & treatment**
- **Pressure relief therapy**
- **Skin care management**

included
**electrical
CONTROL
PUMP**



PVC mattress cover (only for model MAT X3)

Extra strong PVC
with 16 cells



Electrical controlpump
to adjust the firmness
of the mattress

MAT-X3

Images for illustration purposes only. Certain features may be optional.




Opening up new horizons in mobility

Technical data **MAT X1 X2 X3** mattress

Control pump		Mattress X1		Mattress X2		Mattress X3	
Power supply	220V / 50 Hz, 8W	Color	Beige	Color	Blue	Color	Black
Air output	6 liters / minute	Material	PVC	Material	PVC	Material	PVC
Pressure range	25-110 mmHg	Max. capacity	125 kg	Max. capacity	150 kg	Max. capacity	130 kg
Cycle time	5 min. - 10 min.	Size	200 x 90 x 6,5 cm	Size	200 x 90 x 9 cm	Size	200 x 90 x 11,4 cm
Fuse rating	1,5A / 250V	Weight	2 kg	Weight	2 kg	Weight	6,5 kg
Weight	1,2 kg					Coversheet	included

Measurements in cm, weight in kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible.




MAT-X1

PVC bubble mattress



Mattress on bed (illustration purposes only)



MAT-X2

PVC comfort mattress



Mattress on bed (illustration purposes only)



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be



VERMEIREN

Opening up new horizons in mobility

Contact your local agent for more information.

