

Goliat

XXL ROLLATOR



Width	650 mm
Length (unfolded)	750 mm
Seat width	550 mm
Seat height	490 - 560 mm
Seat depth	410 mm
Height handles	830 - 1020 mm
Brakes	multifunctional
Max. inclination	15°
Weight	10,7 kg
Max. users weight	200 kg
Max. capacity basket	5 kg
Wheels	8"
Colours	red / silver

Measurement tolerances: ± 1.5 mm / 1,5kg.
We reserve the right to introduce technical changes.



Model shown only for illustration purposes. Certain features may be optional.

N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: 00 32 (0)3 620 20 20 - Fax: 00 32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com

Opening up new horizons in mobility