

# Bobby EVO

## Lightweight foldable transport wheelchair

**The Bobby EVO is your ideal compact, foldable transfer chair.**  
Light and easy to fold in seconds. Complete with padded, swing away armtubes and a foldable backrest making it easy to stow your chair away.  
Black front and rearwheels and removable and collapsible footrests.



**VERMEIREN**  
we care for you

## Compact and lightweight

The **Bobby EVO** is your ideal compact, foldable transferchair. Light and easy to fold in seconds.

- Swing away armtubes
- 2 manual brakes
- Brakes for the accompanying person
- Padded armtubes
- Angle adjustable footplates
- Adjustable ankle straps

Either on holiday, at the airport or in a hospital, a transfer from point A to point B is easy with the **Bobby EVO**.



## Technical data

■ Total width	490-550 mm
■ Total height	900 mm
■ Total length	1010 mm
■ Folded width	260 mm
■ Seat width	420-480 mm
■ Height armrests	200 mm

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **Bobby EVO** complies with the highest norms, laid down by certified test institutes in most European countries.



■ Seat height	490 mm
■ Height backrest	440 mm
■ Seat depth	420 mm
■ Weight	9,6-10 kg
■ Maximum user weight	115 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer