

708D

Timeless classic



Amongst all the manual wheelchairs available, the 708D is a timeless piece of engineering that has been in demand ever since its first conception. The armrests are foldable and removable and equipped with soft, extra long and comfortable padding. The fabric of the seat and back, together with aluminium leg rests, will provide optimal comfort while using this chair. And as a modular wheelchair, a wide range of options can be offered by **VERMEIREN** to adjust the 708D to your personal taste.



Images for illustration purposes only. Certain features may be optional.



VERMEIREN
we care for you

Standard but complete

Standard wheelchair with foldable seat and backrest. The armrests can be swiveled and removed and are equipped with soft arm pads that are extra long. Complete with padded seat and back as well as aluminium footplates, ensuring a high level of comfort when using the 708D.

Alternate versions



The **T30 version** has small rear wheels that are suitable for usage indoor and outdoor. The wheelchair is intended for users who need to be pushed.

HEM2. With the double grip on one side of the wheelchair, the user can drive and steer his wheelchair with just one hand. The system can be replaced at any time on the desired driving side.



708 D Model S is equipped with an adapter block to lower the seating height for shorter people.

Technical data

| | |
|-----------------|--------------------------------------|
| Total width | 560 / 590 / 610 / 630 / 650 / 670 mm |
| Folded width | 280 mm |
| Seat width | 390 / 420 / 440 / 460 / 480 / 500 mm |
| Height armrests | 225 mm |
| Seat height | 490 mm |
| Height backrest | 420 mm |

Colors



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **708D** complies with the highest norms, laid down by certified test institutes in most European countries.



| | |
|-----------------|---|
| Total height | 910 mm |
| Seat depth | 430 mm |
| Total length | 1109 mm |
| Weight | 18 kg |
| Users weight | 120 kg |
| Additional info |   |

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you